

Shakers Incubators Centrifuges Balances Rotators Vortexers DryBaths Rockers and MORE!



## INCU-SHAKER™ WITH SHAKING INCUBATORS

**MAGiC CLAMP™**  
PATENTED

- Exceptional temperature uniformity
- Large flask capacity
- MAGiC Clamp accessories for instant-magnetic clamp exchange
- Available with CO<sub>2</sub> control or cooling
- NEW! Adjustable clamp for 125, 250 & 500ml flasks

Item: H1000-M-E

**RM 12,500** (Excluded clamps)

INCU-SHAKER Mini

Speed Range: 30 to 300rpm  
Temp Range: Amb. +5 to 70°C  
Platform Dims: 9.5 x 11.5 in.  
Universal Platform: H1000-MR  
Flask Capacity: 2 x 1L (H1000-MR-1000)  
5 x 500ml (H1000-MR-500)  
8 x 250ml (H1000-MR-250)  
12 x 125ml (H1000-MR-125)



Item: H1010-E

**RM 25,500** (Excluded clamps)

INCU-SHAKER 10L

Speed Range: 30 to 300rpm  
Temp Range: Amb. +5 to 60°C  
Platform Dims: 17.8 x 17.8 in.  
Universal Platform: H1010-MR  
Flask Capacity: 5 x 2L (H1000-MR-2000)  
9 x 1L (H1000-MR-1000)  
16 x 500ml (H1000-MR-500)  
25 x 250ml (H1000-MR-250)  
36 x 125ml (H1000-MR-125)



Item: H1012-E

**RM 36,500** (Excluded clamps)

INCU-SHAKER 10LR (with cooling)

Speed Range: 30 to 300rpm  
Temp Range: Amb. -15 to 60°C  
Platform Dims: 17.8 x 17.8 in.  
Universal Platform: H1010-MR  
Flask Capacity: 5 x 2L (H1000-MR-2000)  
9 x 1L (H1000-MR-1000)  
16 x 500ml (H1000-MR-500)  
25 x 250ml (H1000-MR-250)  
36 x 125ml (H1000-MR-125)



Item: H3501-E

**RM 38,500** (Excluded clamps)

INCU-SHAKER Mini CO<sub>2</sub>

Speed Range: 30 to 300rpm  
Temp Range: Amb. +5 to 60°C  
CO<sub>2</sub> Range: 0 to 20%  
Platform Dims: 9.5 x 11.5 in.  
Universal Platform: H1000-MR  
Flask Capacity: 5 x 500ml (H1000-MR-500)  
8 x 250ml (H1000-MR-250)  
12 x 125ml (H1000-MR-125)



Limited Time Offer\*  
**FREE!** RM 2,100 Value  
myTEMP Mini

Item: H2200-H-E  
(Ambient +5 to 60°C)



\*With purchase of any Incu-Shaker  
Special offers are valid through November 15th 2017.



**MC-12™**  
High Speed Microcentrifuge

- Less than 9 inches wide
- LCD with rcf conversion
- Stops in just 10 seconds
- High speed, up to 15,500rpm

**RM 5,300**

Max. Speed: 15,500rpm  
Max. RCF: 16,100xg  
Tilt Angle: 12 x 1.5ml

C1612-E MC-12 High Speed Microcentrifuge with 12 x 1.5ml rotor RM 5,300

Limited Time Offer\*  
**FREE!**

RM 1,600 Value



**myFUGE 12**  
Mini Centrifuge

\*With Purchase of C1612 or B3D5000

**BENCHWAVER™**  
Undulating 3-D Shaker

- Undulating 3D motion
- Adjustable speed and tilt
- Ideal for staining, washing and general mixing



**RM 4,900**

Speed: 5-105rpm  
Platform: 13.5 x 13 in.  
Tilt Angle: 0 to ±10°

Optional Stacking Platform

B3D5000-E BenchWaver™ Undulating Shaker RM 4,900  
B3D5000-STK Optional stacking platform RM 500  
B3D5000-DIMP Optional dimpled mat RM 400

